

IL ADT integration setup

Go to API settings on ADT tools /ADT/tools/api/settings and follow the below steps

ADT v3.4(Beta-1.0) Login

RecoverBackupMigrationSystem UpdateDB ManagerHelpInitial SetupAPI Settings

API Settings

1

Api Status

2

Api Patients Module

Api Dispense Module

Api Appointments Module

3

Api Adt Url192.168.0.223/ADT/tools/api

Api Adt Port6671

Api Il Ip192.168.0.9

Api Il Port9720

4

Api Logging

Save Settings

Government of Kenya © 2018. All Rights Reserved

1. Turn on API

Api Status ☒

2. Turn on ADT modules you wish to communicate with IL

Api Patients Module ☒

Api Dispense Module ☒

Api Appointments Module ☒

3. Specify ADT and IL urls and ports

Api Adt Url	192.168.0.223/ADT/tools/api
Api Adt Port	6671
Api Il Ip	192.168.0.9
Api Il Port	9720

4. Turn on API logging

Api Logging ☒

To see whether messages are being generated go to /ADT/tools/IL-api.log

Logs are written in the below format

YYYY-MM-DD H:M:S {message type}: message

e.g 2018-05-23 12:35:14 PATIENT ORU^R01 :
{ "MESSAGE_HEADER": { "SENDING_APPLICATION": "ADT", "SENDING_FACILITY": "13050", "RECEIVING_APPLICATION": "IL", "RECEIVING_FACILITY": "13050", "MESSAGE_DATETIME": "20180523123514", "SECURITY": "", "MESSAGE_TYPE": "ORU^R01", "PROCESSING_ID": "P" }, "PATIENT_IDENTIFICATION": { "EXTERNAL_PATIENT_ID": { "ID": "", "IDENTIFIER_TYPE": "", "ASSIGNING_AUTHORITY": "" }, "INTERNAL_PATIENT_ID": [{ "ID": "1841", "IDENTIFIER_TYPE": "SOURCE_SYSTEM_ID", "ASSIGNING_AUTHORITY": "ADT" }, { "ID": "13050-PEP428", "IDENTIFIER_TYPE": "CCC_NUMBER", "ASSIGNING_AUTHORITY": "CCC" }], "PATIENT_NAME": { "FIRST_NAME": "WINNIE", "MIDDLE_NAME": "SILAS", "LAST_NAME": "" }, "OBSERVATION_RESULT": [{ "SET_ID": 1, "OBSERVATION_IDENTIFIER": "START_HEIGHT", "CODING_SYSTEM": 1, "VALUE_TYPE": "NM", "OBSERVATION_VALUE": "null", "UNITS": "CM", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }, { "SET_ID": "2", "OBSERVATION_IDENTIFIER": "START_WEIGHT", "CODING_SYSTEM": "", "VALUE_TYPE": "NM", "OBSERVATION_VALUE": "null", "UNITS": "KG", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }, { "SET_ID": "3", "OBSERVATION_IDENTIFIER": "IS_PREGNANT", "CODING_SYSTEM": "", "VALUE_TYPE": "CE", "OBSERVATION_VALUE": "NO", "UNITS": "YES\\NO", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }, { "SET_ID": "4", "OBSERVATION_IDENTIFIER": "PRENANT_EDD", "CODING_SYSTEM": "", "VALUE_TYPE": "D", "OBSERVATION_VALUE": "20170713110000", "UNITS": "DATE", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }, { "SET_ID": "5", "OBSERVATION_IDENTIFIER": "CURRENT_REGIMEN", "CODING_SYSTEM": "NASCOP_CODES", "VALUE_TYPE": "CE", "OBSERVATION_VALUE": "38", "UNITS": "", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }, { "SET_ID": "6", "OBSERVATION_IDENTIFIER": "IS_SMOKER", "CODING_SYSTEM": "", "VALUE_TYPE": "CE", "OBSERVATION_VALUE": "NO", "UNITS": "YES\\NO", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }, { "SET_ID": "6", "OBSERVATION_IDENTIFIER": "IS_ALCOHOLIC", "CODING_SYSTEM": "", "VALUE_TYPE": "CE", "OBSERVATION_VALUE": "NO", "UNITS": "YES\\NO", "OBSERVATION_RESULT_STATUS": "F", "OBSERVATION_DATETIME": "20130115120000", "ABNORMAL_FLAGS": "N" }] }

5. Save settings

Click on Save settings button to save your changes.

Save Settings

Once you have saved settings your page should look like this

ADT v3.4(Beta-1.0) Login

RecoverBackupMigration ▼System UpdateDB ManagerHelpInitial SetupAPI Settings

API Settings

Api Status

☒

Api Patients Module

☒

Api Dispense Module

☒

Api Appointments Module

☒

Api Adt Url

192.168.0.223/ADT/tools/api

Api Adt Port

6671

Api Ii Ip

192.168.0.9

Api Ii Port

9720

Api Logging

☒

Save Settings